



Philips
Shaving heads

HQ4

Renew for better results

For maximum performance replace your shaving heads every 2 years

Microgroove

- Microgroove for a closer shave

Lift & Cut shaving system

- Dual blade system lifts hair to cut comfortably

15 razor-sharp blades

- For a fast shave

PHILIPS
sense and simplicity

Specifications

Shaving heads

- **Shaving heads per packaging:** 3
- **Fits product types:** HQ4401, HQ4405, HQ4407, HQ4411, HQ443, HQ4445, HQ4800, HQ4806, HQ4807, HQ4810, HQ4819, HQ482, HQ4826, HQ4830, HQ484, HQ4846, HQ485, HQ4850, HQ4856, HQ486, HQ4870, HQ4890, HQ5812, HS 190, HS766, HS775, HS890, HS920, HS930, HS970, HS990



Issue date 2009-02-10

Version: 3.1.1

12 NC: 0000 000 00000

© 2009 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips
Electronics N.V. or their respective owners.

www.philips.com