

Philips
Digital thermometer set

SCH530/86



Digital body thermometer

Measure within seconds

Dimpled shield on soother

Silicon teat

Detachable electronics

Easy and reliable temperature measurement

Digital thermometer soother for easy and accurate measurement.

Fast and reliable measurement

- Digital thermometer for easy and accurate measurement
- Quick and reliable temperature measurement within seconds

Convenience for you and your baby

- Dimpled shield prevents skin irritation
- Designed especially for newborn babies
- Includes a dummy pacifier in the size for a newborn baby

Baby's safety

- Allows boiling or sterilising without damaging electronics
- Teat suitable for sterilising or boiling

PHILIPS
sense and simplicity

Specifications

Technical specifications

- **Operating temperature range:** 32–44 °C
- **Accuracy:** 36.5–39.0C:±/0.1C °C
- **Measurement time:** 10 s

Convenience

- Automatic shut-off
- Ready indicator
- Flexible tip
- Suitable for newborns

Accessories

- Battery
- User manual

Power

- **Battery type:** LR41
- **Number of batteries:** 2
- **Operating time on battery:** 200 hour(s)
- Removable/replaceable

Product highlights

Digital body thermometer

Fast, reliable and accurate digital thermometer for measuring the baby's temperature.

Measure within seconds

Never too hot, never too cold — the "Just Right" digital thermometer quickly tells you if the milk or food is safe to eat or drink. It takes the guessing out of warming bottles or jars and is great for checking if food is still warm enough after taking a short break in feeding.

Dimpled shield on soother

Dimples make sure air flows behind the teat to protect your baby's sensitive skin from irritation.

Silicon teat

The soother's specially shaped silicon teat is designed to fit your baby's mouth and will not inhibit natural development. Use it with confidence to pacify your baby.

Detachable electronics

By detaching the sensitive temperature measuring electronics from the soother, you can safely sterilise or boil the part your baby comes into contact with.

Sterilise teat

Thanks to the special material, you can safely sterilise the soother's teat in a microwave oven or electrical steriliser or boil it without harming the soother.

