

# PHILIPS

## Saeco

### Espresso machine descaler

for Saeco Espresso machines

1 descaling cycle

reduces limescale

CA6700/99



## Perfect protection for your espresso machine

### Special Saeco espresso machine descaler

Decalcifying your espresso machine regularly is essential to make sure it keeps performing at its best. This special espresso machine descaler removes scale and prevents corrosion, protecting your appliance and extending its lifetime.

#### Exclusive formula

- For thorough yet safe and gentle decalcification
- Refined over 25 years

#### Effective protection of your espresso appliance

- Perfect decalcification of water circuits
- Protects your system against limescale build-up

#### Longer lifetime

- Regular decalcification prolongs lifetime and enhances taste

#### Quick and easy decalcification

- The best formula for easy decalcification

#### Sensational coffee quality

- Maintains the coffee taste over time

# Highlights


## Exclusive formula

The exclusive formula of the Philips Saeco espresso machine decalcifier ensures thorough decalcification without harming the delicate parts inside your appliance.

## Refined over 25 years

Tested and approved by Saeco, the inventor of the full-automatic espresso machine.

## Special Espresso Decalcifier



The Philips Saeco decalcifier perfectly cleans all the water circuits in your espresso machine.

## Reduces limescale build-up

Limescale is a natural part of water used for the machine's operation. This special decalcifier protects your appliance from the

build-up of limescale which affects performance and taste. It is highly effective, safe and easy to use. For effective decalcification descale when the machine requires it or after 250 cups according to water hardness.

## Better taste and aroma



Regular maintenance ensures the best taste and aroma from your Saeco espresso machine.

## Quick and easy decalcification

The special Saeco decalcifier for espresso machines comes with clear instructions for use, enabling quick, easy and safe decalcification.

# Specifications

## Country of origin

Made in: Germany

## Technical specifications

Applications: 1 descaling cycle

Contents: 1 bottle of 250 ml

Ingredients: 50% citric acid, 25% lactic acid

## Weight and dimensions

Product dimensions (L x D x H): 66 x 40 x 169 mm

Product weight: 0,25 kg

