

**PHILIPS**

**AVENT**

Silicone diaphragm  
for breast pumps

CP9823



To bring vacuum in the breast pump.

Find matching products on the specifications tab

This silicone diaphragm creates vacuum in your manual breast pump.

**Ease of using Philips spare parts**

- Easy renew your product with original Philips parts

## Highlights

### Easy renew of product parts

From time to time your product needs a facelift but it was never so easy like with Philips consumer replacement parts to renew your product. This all with guaranteed Philips quality.

## Specifications

### Replaceable part

Fits product types :: SCF330

